

## Features

- 350 lb. weight capacity
- 20 miles estimated range
- 12 degree maximum incline
- 0-4 mph variable speed
- 31" turning radius
- Adjustable position, programmable Dynamic controller with common connectors
- Removable, flip-up armrests with adjustable height, width and angle
- Flip-up footplate with adjustable height and angle
- 18", 20" & 22" seat widths
- Semi-reclining, foldable Captain's Seats with headrests and 180-degree swivel
- Batteries included

## Specifications

<b>Max. Weight Capacity</b>	350 lbs.
<b>Maximum Speed</b>	4 mph
<b>Estimated Range</b>	20 miles
<b>Maximum Grade/Incline</b>	12°
<b>Turning Radius</b>	31"
<b>Ground Clearance</b>	3"
<b>Overall Dimensions</b>	
Length	43"
Width	24"
Height	42"
<b>Power</b>	
Motor	Two 24V x 200W x 4600rpm
Controller	Dynamic A50A
Batteries	Two DC12V x 36AH
Battery Charger	5 Amp, Offboard
<b>Brakes</b>	Electromagnetic
<b>Freewheel Mode</b>	Yes
<b>Tires</b>	
Front Casters (Foam-Filled)	8" x 2"
Rear Drive Wheels (Pneumatic)	10" x 3"
<b>Weights</b>	
Without Batteries	143 lbs.
With Batteries	191 lbs.




## Accessories



## Available Colors



## Seating Options

<b>Renegade 18</b>	
	Width 18"
	Depth 18"
	Back 29"
	Armrest to Seat 7"
<b>Renegade 20</b>	
	Width 20"
	Depth 20"
	Back 29"
	Armrest to Seat 7"
<b>Renegade 22</b>	
	Width 22"
	Depth 20"
	Back 29"
	Armrest to Seat 7"

**ActiveCare**  
Medical

866.621.4933 | [activecaremed.com](http://activecaremed.com)