

## Features



Spitfire I410

- Available in 3-wheel (I310) and 4-wheel (I410) versions
- Highest-in-class 300 lb. weight capacity
- Disassembles without tools into 5 pieces, the heaviest weighing only 28 pounds
- 9 mile maximum range
- Adjustable speed up to a maximum of 5 mph
- Adjustable angle tiller with fingertip throttle control
- Quick connect battery pack
- Anti-tip wheels for greater safety
- Flip-up armrests with adjustable width and angle
- Anti-tip wheels for greater safety
- Comfortable, 17" wide swivel seat with fold-down, padded backrest
- Standard metal wire basket
- Flat-free, non-marking tires
- Batteries included



Spitfire I310



Spitfire scooters quickly and easily disassemble into 5 pieces

## Color Options



## Accessories

- Oval Bag
- Carry All Bag
- Dust Cover

## Specifications

<b>Max Weight Capacity</b>	<b>300 lbs.</b>	<b>Brakes</b>	Electromagnetic
<b>Max Speed</b>	5 mph	<b>Freewheel Mode</b>	Yes
<b>Max Range</b>	9 miles	<b>Tires (Flat-free)</b>	
<b>Max Grade/Incline</b>	12°	Front & Rear	8" x 2"
<b>Turning Radius</b>	43" (I310)   47" (I410)	<b>Weights</b>	
<b>Ground Clearance</b>	2"	Without Batteries	66 lbs. (I310) 77 lbs. (I410)
<b>Dimensions</b>		With Batteries	88 lbs. (I310) 99 lbs. (I410)
Overall Length	42" (I310)   43" (I410)	Heaviest Piece Wt.	28 lbs. (rear section)
Overall Width	19"		
Floor to Seat Height	19"-22"		
<b>Seating Dimensions</b>			
Seat Size	17"(w) x 16"(d)		
Back Height	14"		
Armrest to Seat Height	7"		
<b>Power</b>			
Motor	24V x 270W x 4700rpm		
Controller	Dynamic R50A		
Batteries	Two DC12V x 12AH		
Battery Charger	2 Amp, Offboard		

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

# ActiveCare

Medical

2 Harbison Way  
Columbia, SC 29212  
866.621.4933  
activecaremed.com

Information contained herein is correct at the time of publication. ActiveCare reserves the right to alter specifications without notice.