


ActiveCare
drive™ Mobility
COBALT X23

NEW!



**FULL SIZE
 PERFORMANCE
 IN A LIGHTWEIGHT
 POWER CHAIR**



Midnight Blue Metallic Red

**PWC GROUP 2
 HCPCS
 CODE:
 K0823**

Cobalt X23

Standard Power Wheelchair

Available
1 Pack

COBALTX23BL16FS
 COBALTX23RD16FS

16", Folding Seat, Midnight Blue, 1/cs
 16", Folding Seat, Metallic Red, 1/cs

FEATURES:

- Full-size performance in a lightweight, easily transportable power wheelchair.
- Seats have four height adjustments, from 20" to 23".
- Comfortably padded, flip-up armrests allow for easy transfers and access to user. Armrests are also removable and height and width adjustable for positioning, support and adaptability to the size of the user.
- Adjustable length controller mount can be affixed to either the left or right armrest and can be fitted with an optional swing-away arm for additional flexibility.
- Standard footplate can be flipped up out of the user's way and adjusted in both height and angle.
- Chair can be fitted with optional elevating legrests allowing the user to travel comfortably with their legs raised—a necessity for those with some medical conditions, such as edema.
- Standard positioning belt for added user security.
- Rear anti-tip wheels for added safety.
- Armrest-mounted reflectors improve the chair's visibility.
- Extra strong steel frame for long-lasting durability.
- Programmable controller with common connectors is easy to service and setup and can be personalized for individuals with special needs.
- Flat-free, non-marking tires offer worry free travel.
- Motor assemblies offer easy "freewheel" operation, allowing the chair to be pushed, which can be helpful for transport.
- 250 lbs. Weight Capacity

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.