

Heartway/ Passport(Lithium)

Model: S19/Z-16 w/ Lithium

Item: 14282028

S19 Passport Folding Scooter (Bronze) -Lithium Battery

Description:

S19 is an ultra light weight folding scooter. It weighs only 54 lbs without battery and so is an ultimate product for traveling carry on mobility equipment. It travels a range of 12.4 miles. It is equipped with a Lithium Battery providing it with ultra light weight and bigger range of drive. It folds up as easy as just one step. Leg room of a 3 wheeled scooter combined with the stability of 4 wheels.



MSRP: USD \$2,995

Product Video: <http://youtu.be/iIo9D0R1bls>

Heartway reserves the rights to alter any content stated above without prior notice.

ITEM	SPECIFICATION
WEIGHT CAPACITY	250 lbs (115kgs)
DRIVE WHEEL	8" x2" (200 x 50mm)
FRONT CASTER	7"x 1.6" (180x 40mm)
REAR CASTER	1.5"x0.75" (38 x 19mm)
MAX SPEED	4 MPH (6.2 KPH)
BATTERY RANGE	12.4 miles (20 km)
CHARGER TYPE	2.5A Off Brd Li-on 110~240V
CONTROLLER TYPE	S-Drive 45A
MOTOR TYPE	4-Pole Motor 270 W
DIMENSIONS:	Unfolded: 37.2"H x19"Wx 36.6"D Folded: 17.7" H x 19" W x 29.5" D
WEIGHT:W/O BATT	54 lbs (24.5kgs)
BATTERY PACK	Lithium 11.4 Ah – 5.5 lbs (2.5 kgs)
SUSPENSION	Front: INDEGO Suspension
TURNING RADIUS	32.2" (817 mm)
SEAT STYLE	Foldable Fish-On Seat
SEAT WIDTH	16" (406 mm)
SEAT HEIGHT	12.6" (320 mm)
SEAT DEPTH	13" (330 mm)
BACK HEIGHT	14.2" (360 mm)
WHEEL BASE	27.4" (696 mm)
GROUND CLEARANCE	2.4" (60.9 mm)
LEG ROOM	16.5" (419 mm)

Standard Package (included):

Battery: 11.5 Ah Lithium Battery x 1

Charger: Off Board -2.5Amp Lithium 110~240V

Seat: 16 in Fish-on seat

Tires: Flat Free Tires

Accessory: Charger, Battery Box, User Manual.

Available Color: Bronze

Free Shipping within 48 U.S. continental states.

Options:

Seat Size: n/a

Tires: n/a

Battery: 14.5 A Lithium call for pricing

Armrest Assembly (+\$180)

Color: Blue

One Year In-Home Limited Labor Svc

HWL Program (+299)

White Glove Delivery (+300)

Accessories:

- Rear Large Basket n/a

- Cane Holder (+\$80)

- Canopy & Rear Bag n/a

- Weather Proof Cover-S (+\$90)

- Oxygen Tank Holder n/a

- Monster Rear Bag n/a

- Locking Trunk Box w/ Key n/a

- Extra Long Seatbelt 72inch n/a

-European Power Adaptor (+\$30)

-Additional Li- Battery Pack (+\$300)

Order: online.heartwayusa.com

Heartway/ Passport (SLA)

Model: S19/Z-16 w/ SLA

Item: 14282028 SLA

S19 Passport Folding Scooter (Bronze) -SLA (AGM) Battery

Description:

S19 is an ultra light weight folding scooter. It weighs only 54 lbs without battery and so is an ultimate product for traveling carry on mobility equipment. With SLA battery, it travels a range of 8 miles. It is equipped with two 12v 12Ah SLA Batteries for a more economy choice. It folds up as easy as just one step. Leg room of a 3 wheeled scooter combined with the stability of 4 wheels.



MSRP: USD \$2,795

Product Video: <http://youtu.be/iIo9D0R1bIs>

ITEM	SPECIFICATION
WEIGHT CAPACITY	250 lbs (115kgs)
DRIVE WHEEL	8" x2" (200 x 50mm)
FRONT CASTER	7"x 1.6" (180x 40mm)
REAR CASTER	1.5"x0.75" (38 x 19mm)
MAX SPEED	4 MPH (6.2 KPH)
BATTERY RANGE	8 miles (12.8 km)
CHARGER TYPE	2A Off Brd (SLA) 110~240V
CONTROLLER TYPE	S-Drive 45A
MOTOR TYPE	4-Pole Motor 270 W
DIMENSIONS:	Unfolded: 37.2"H x19"Wx 36.6"D Folded: 17.7" H x 19" W x 29.5" D
WEIGHT:W/O BATT	54 lbs (24.5kgs)
BATTERY PACK	TWO 12V 12Ah – 17 lbs (7.7 kgs)
SUSPENSION	Front: INDEGO Suspension
TURNING RADIUS	32.2" (817 mm)
HEIGHT	34" (860 mm)
SEAT WIDTH	16" (406 mm)
SEAT HEIGHT	12.6" (320 mm)
SEAT DEPTH	13" (330 mm)
BACK HEIGHT	14.2" (360 mm)
WHEEL BASE	27.4" (696 mm)
GROUND CLEARANCE	2.4" (60.9 mm)
LEG ROOM	16.5" (419 mm)

Standard Package (included):

Battery: 12Ah Lead-Acid (SLA) Battery x 2pcs
Charger: Off Board -2Amp Charger 110~240V
Seat: 16 in Fish-on seat
Tires: Flat Free Tires
Accessory: Charger, Battery Box, User Manual.
Available Color: Bronze

Free Shipping within 48 U.S. continental states.

Options:

Seat Size: n/a
Tires: n/a
Battery: 14.5 A Lithium call for pricing
Armrest Assembly (+\$180)
Color: Blue

One Year In-Home Limited Labor Svc
HWL Program (+299)

White Glove Delivery (+300)

Accessories:

- Rear Large Basket n/a
- Cane Holder (+\$80)
- Canopy & Rear Bag n/a
- Weather Proof Cover-S (+\$90)
- Oxygen Tank Holder n/a
- Monster Rear Bag n/a
- Locking Trunk Box w/ Key n/a
- Extra Long Seatbelt 72inch n/a
- European Power Adaptor (+\$30)
- Additional Li- Battery Pack (+\$300)

Order: online.heartwayusa.com