

Jazzy[®]

SELECT ELITE

The Jazzy Select[®] Elite delivers a potent blend of power, performance and style. In-line, front-wheel drive technology gives the Elite excellent stability and maneuverability for solid performance indoors and out.



Jazzy Select[®] Elite shown with new gray vinyl comfort seat, PG GC 2 controller and foot platform

800-800-8586 (US) • 888-570-1113 (Canada)
Visit us on the web: www.pridemobility.com

Pride[®]
Mobility Products Corp.

SPECIFICATIONS

Weight capacity	300 lbs.
Maximum speed	Up to 4 mph
Ground clearance	1.5" at anti-tips
Turning radius	24.75"
Overall length	34.75" without front riggings
Overall width	22.75"
Wheels	Drive: 9"
Front anti-tips	3" solid
Rear casters	6" solid
Suspension	Limited
Drive train	Two-motor, in-line, front-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC 2 controller
Specialty controls	N/A
Battery charger	Off-board, 3 amp
Per-charge range	Up to 15 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight	24.5 lbs. each
Base weight	68.5 lbs.
Standard seat weight	43.5 lbs.
Max. seat dimensions	20" x 20"
Warranty	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors

COLOR CHOICES



Jazzy Red



Jazzy Blue

MODELS

Jazzy Select® Elite-2S-C
Jazzy Select® Elite-2S-SS

FEATURES

- In-line, front-wheel drive technology
- 40 amp, PG GC2 controller
- Easy access to batteries from rear
- Simple main frame design for easy serviceability
- New comfort high-back seat with headrest

OPTIONS

- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

RECOMMENDED LIFTS

Silver Star®
Lifts & Ramps

- Backpacker® Plus
- Full Platform
- Backpacker® MV



800-800-8586 (US) ▪ 888-570-1113 (Canada)
Visit us on the web: www.pridemobility.com

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

PRIDE®
PERFORMANCE MOBILITY

— LIVE YOUR BEST™ —