

# Luggie Standard



## Compare Models

Features	Classic	Standard	Elite
Weight Capacity	250 Lbs	250 Lbs	320 Lbs
Range*	10 Miles	11 Miles	14 Miles
Battery	Sealed Lead Acid	8.5 AH Lithium	10.5 AH Lithium
Battery Weight	20 Lbs	5 Lbs	5 Lbs
Battery Charging Time	4-6 Hrs	2 Hrs	2 Hrs
Seat Width	14.5 inches	14.5 inches	16 inches
Max Speed	4 MPH	4 MPH	4 MPH
Battery Indicator	Meter	Meter	LED
Speed Control Dial	Yes	Yes	Yes
Extendable Handlebars	No	No	Yes
Overall Length (unfolded)	39"	39"	39"
Overall Length (folded)	26"	26"	26"
Overall Width	18"	18"	18"
2 Seat Height Options (20"&22")	Yes	Yes	Yes
Turning Radius	35.5"	35.5"	35.5"
Ground Clearance	2.5"	2.5"	2.5"
Adjustable Steering Heights	Yes	Yes	Yes
Auto Retractable Bumper	No	No	Yes
Total Weight (without battery)	51.7 Lbs	51.7 Lbs	51.7 Lbs
Aluminum Alloy Frame	Yes	Yes	Yes
Rear wheel drive	Yes	Yes	Yes
24 Volt DC magnetic motor	Yes	Yes	Yes
Intelligent Braking System	Yes	Yes	Yes
Airline Safe	Yes	Yes	Yes