

## P320 Junior

	<p>Technical Specs</p> <p><b>Max Speed:</b> 4.5mph  <b>Range:</b> up to 12.5 miles  <b>Turning Radius:</b> 23"  <b>Ground Clearance:</b> 2"  <b>Motor:</b> DC 24V, 120W  <b>Controller:</b> Dynamic Shark 60A/PG VR2 40A  <b>Battery:</b> 12V/22A x 2pcs  <b>Charger:</b> 3A off-board  <b>Gradient:</b> 8°  <b>Casters:</b> 7 x 1.75" solid tires  <b>Drive Wheels:</b> 9" foam filled tires  <b>Brakes:</b> intelligent, regenerative, electromagnetic</p>
<p>MSRP</p> <p>US: P320-1S C =\$4,499.00  P320-1S S  P320-2S C =\$4,499.00  P320-2S S</p>	
<p>Medicare Code</p> <p>K0816 P320-1S C  K0815 P320-1S S  K0823 P320-2S C  K0822 P320-2S S</p>	
<p>Certification</p> 	
<p>Features</p> <ul style="list-style-type: none"> <li>• Lightweight and easy to disassemble for handling, storage and transportation</li> <li>• Ideal as personal secondary portable power wheelchair</li> <li>• All standard power wheelchair features</li> <li>• Programmable controller with a wide range of selections to suit the needs and demands of different users</li> <li>• Non-programmable controller also available</li> <li>• Excellent maneuverability</li> <li>• Off-board charger</li> <li>• Colors Available: Red, Blue</li> <li>• Options: Elevating Leg Rests, Oxygen Tank Holder, Cane &amp; Crutch Holder, Walker Holder, Foldable Rear Basket</li> </ul>	

## Dimensions

Seat Width	18
Seat Depth	16
Seat Height	16-18
Armrest Height	n/a
Backrest Height	17
Min. Footrest Extension	n/a
Max. Footrest Extension	n/a
Overall Length	33
Overall Height	35
Overall Width	24
Folded Width	n/a
Total Weight	135 lbs.
Capacity	300 lbs.