



P326A Vision Sport - Mid Wheel Drive Power Wheelchair

	<p>Technical Specs</p> <p>Max Speed: 5mph Range: up to 18 miles Turning Radius: 20" Ground Clearance: 2" Motor: DC 24V, 160W Controller: Dynamic Shark 50A Battery: 12V/U1 x 2pcs Charger: 4A off-board Gradient: 6° Casters: 6" front and rear solid tires Drive Wheels: 10" foam filled tires Brakes: intelligent, regenerative, electromagnetic</p>
<p>MSRP</p> <p>US: \$6,399.00</p>	
<p>Medicare Code</p> <p>P326-2S S (K0822) P326-2S C (K0823)</p>	
<p>Certification</p> 	
<p>Features</p> <ul style="list-style-type: none"> • Full suspension • No rear tipping on footplate standing (up to 250 lbs.) • 300 lb weight capacity • Dual in-line motors for enhanced efficiency, torque, range and performance • Swing back joystick mount is standard • Height and width adjustable arms • Semi-reclining back • Optional Accessories: Elevating legrests, Oxygen tank holder, Cane & Crutch holder, Walker holder, Foldable rear basket, Canopy and Variety of Color Fenders 	

Dimensions

Seat Width	20
Seat Depth	18

Seat Height	21.5-23.5
Armrest Height	n/a
Backrest Height	20
Min. Footrest Extension	n/a
Max. Footrest Extension	n/a
Overall Length	34
Overall Height	48.5
Overall Width	24
Folded Width	n/a
Total Weight	184 lbs.
Capacity	300 lbs.