

PaceSaver®
Set Your Own Pace!

Scout™ BOSS

6.75 & 6

Midi-Drive™ Power Wheelchairs

Attractive BOSS power bases put you in control.

BOSS 6.75 & Boss 6 Excel with ...

Style & choice - 3 rich "Mirror-Glaze™" Through & Through automotive colors, 4 different seat colors, including 2 luxurious soft fabrics.

Easy chair access - just step on the footrest.

All day comfort - Orthopedic contoured Com-For-Back™ seat and arms are width, depth, height & angle adjustable to fit your needs. 50° Recline Admiral Seat. SS is available in Com-For-Back™, Admiral & Rehab Cane Seat.

Better indoor maneuverability and stability - Patented 3/4 midi-drive short 27" turning radius is much easier to drive in your home than comparable rear wheel drive chairs. Boss' 3/4 design is more stable than center drive bariatric chairs that tip easily when reaching or stopping. 8° stability.

All day outdoor power & range - New York made, (1^{3/4} HP Total) high torque motors drive maximum efficiency, inline (no maintenance) sealed gear drive. Advanced, easy to drive, (180/6.75 150/6) amp modular controller saves energy while delivering maximum power for range up to (18/6.75 - 20/6) miles. Boss 6.75/Boss 6 can climb tortuous, 140yd, 7-9° hills 20/24 times with a 675lb./600lb. rider.

Shock absorbing spring suspension for smoother ride while climbing 2 1/2" obstacles, rocks or gravel.

Boss 6.75's / 6's ultra low center of gravity is designed with 4 point seat mounts that ensure much lower seat heights, (16" for a short user - SS), for better control, drivability and easy battery access.

Fast charging anywhere. Small carry-along worldwide charger with new state-of-the-art battery saver technology for longer battery life.

Stability & traction - 3" seat to base balance adjustment optimizes drive for any size or shape.

Transportability, separates into 6 pieces for transport or use chair tie down system for in-vehicle transport.

Special Need Extra wide footplate standard. Swing-away, elevating legrests optional.

Made in America for custom seating & quick parts availability.



Shown with optional 50° Recline Admiral Seat and adjustable headrest extension



675lb. / 600lb. capacity

Meet The "BOSS 6.75 & BOSS 6"

PaceSaver Celebrating 23 Years of legendary performance

PaceSaver®
Set Your Own Pace!

Scout™ BOSS
6.75 & 6
Midi-Drive™ Power Wheelchairs

ALWAYS
AMERICAN MADE



Com-For-Back™ Seat (CC)

Back folds for transport

20° back angle adjustment (16° posterior, 4° anterior)
Back Height 22" or 25" with 5" or 8" posterior space
3" back frame depth adjustment

Deluxe "Com-For-Back™" seat

24" width standard (26" to 34" optional)
21" depth standard (custom sizes available)
3" height adjustment, 5.4° tilt angle adjustment
Top of seat to floor: 18.3" - 21.3" (uncompressed foam)
(Taller heights available)
3" seat on frame position adjustment (front to rear)

3" x 14" Removable arms

Width, height, depth & angle adjustable
3" depth adjustment

22" inside width, 28.5" outside width to 30" inside width
with 36" outside width

Removable Footrest

3" to 5" height adjustable and angle adjustable



1 of 4 Solid Seat Styles (SS)

Specifications subject to change without notice.

Primary Use	Indoor/Outdoor
Speed	4.5 mph
Range (varies with weight)	6.75 (18 miles / 675 lbs.) / 6 (20 miles / 600 lbs.)
Ground Clearance	3"
Turning Radius	27"
Overall Length	48"
Overall Width	28.5"
Overall Height	42" - 44"
Batteries (2 Required)	Group 24 Sealed
Battery Charger	Carry-On, 8 amp, 24 volt constant current
Weight	
Main Frame & Components	147 lbs.
Footplate/Front Anti-Tip System	26 lbs.
Seat and Controller	59 lbs.
Batteries (2)	58.5 lbs. each
TOTAL	349 lbs. includes batteries & flat free tires
Structural Capacity	6.75 (675 lbs.) / 6 (600 lbs.)
Maximum Incline Rating	8°
Controller - Shark Controller	Fully Programmable
Total Amps	6.75 (180 amps) / 6 (150 amps)
Suspension/System	
Main	Coil Spring
Front Anti-Tip System	Fully Adjustable
Tires & Wheels	
Main Drive	12" x 3" Flat Free (Stn) or (Optional Pneumatic -NC)
Front	5" Solid
Rear	9" x 2.5" Flat Free (Stn) or (Optional Pneumatic -NC)

Warranty: Lifetime limited warranty on structural frame; two-year limited on electronics; two year limited on motor and gearbox assemblies. Please see Owners Manual for details.



Available with optional
Stanza™ Power Anti-tips
for easy step on footrest
access and extra ramp stability.

Accessories & Options

- Front Bag
- Folding Basket
- Cup Holder
- Swing Away Controller Mount
- Adjustable Head Rest
- Head Rest Extension
- Oxygen Cylinder Holder
- Elevating Leg Rests
- Swing-Away/Detachable Leg Rests
- Saddle Bag
- Cane Holder
- Rehab (cane) Seat

Body Color Choices

- Ocean Blue
- Candy Apple Red
- Midnight Black

Seat Color Choices - Velvet

- Black
- Gray

Seat Color Choices - Vinyl

- Gray
- Charcoal

PaceSaver® Celebrating 23 years of legendary performance!

PaceSaver® and Leisure-Lift™ are registered trademarks and Scout™, Boss™, Com-For-Back™, Midi-Drive™, and Stanza™ are trademark of Leisure-Lift, Inc. Shark™ is a trademark of Dynamic Controllers Ltd. ©2005 Leisure-Lift, Inc. All rights reserved. Numerous Patents Pending

www.pacesaver.com
leisure-lift@kc.rr.com

Leisure-Lift, Inc.
1800 Merriam Lane
Kansas City, KS 66106
1-800-255-0285
Boss675.indd