

# Solax Mobie



## Specification

|                              |   |
|------------------------------|---|
| Weight Capacity:             | 300 lbs.  |
| Total Weight - w/ Batteries: | 54 lbs.   |
| Total Weight - w/o Batteries | 50lbs.  |
| Freewheel Mode:              | Yes   |
| Maximum Speed:               | 3.75 mph  |
| Estimated Range:             | 13.5 miles                                      |
| Maximum Incline:             | 8 degrees                                       |
| Turning Radius:              | 47"   |
| Ground Clearance:            | 1.6"  |
| Overall Length:              | 38.5"   |
| Overall Width:               | 17.3"   |
| Overall Height:              | 20"   |
| Seat Width:                  | 18"   |
| Seat Depth:                  | 16"   |
| Seat Back:                   | 14"   |
| Armrest to Seat:             | 7"  |
| Folded Dimensions:           | 23.5" H, 17" W, 16" D                           |
| Motor:                       | 120W  |
| Batteries:                   | 24V / 10Ah (Lithium Battery)                    |
| Battery Charger:             | 13.5 miles                                      |
| Brakes:                      | Intelligent regenerative Electromagnetic Brakes |
| Tires:                       | 8" x 2" Flat Free                               |