

Solax Transformer



Features

The Transformer is the only fully electric folding scooter in the USA. There is nothing else like it on the market. This easy to use, ultra portable electric scooter offers the ultimate in convenience. Offering Remote Control operation and manual if required. The Transformer has a 300lb weight capacity and comes standard with a lightweight Lithium battery (13.5 Miles on a single charge).

- Remote Control
- Fully Automatic Folding & Unfolding
- Ultra Portable
- Perfect for Aircraft & Cruise Ships
- 300lb Weight Capacity
- 13.5 Miles on a Single Charge
- Lithium Battery

Product Specifications

| | |
|------------------------------------|---|
| Model: | S3021 Automatic folding scooter |
| Overall Length: | 0.93m (36.9in) |
| Overall Width: | 0.45m (18.1in) |
| Tires (Solid wheels): | Tire Size Front 6 "x 1.5" Tire Size Rear 7 "x 2.36 " |
| Maximum Speed: | 6km/h (3.7mph) |
| Safe Gradient/ Maximum Gradient | 0-15° |
| Range Up To: | 20km (13.7mile) |
| Turning Circle: | 1.4m (55in) |
| Weight (without battery pack): | 24kg (53lbs) |
| Color: | Red, yellow, blue |
| Brake: | Intelligent, regenerative, and |

| | |
|--------------------------------|-----------------------------|
| | electromagnetic |
| Drive system: | S-drive 45A |
| Handlebar: | PU foam |
| Maximum weight capacity: | 136kg (300lbs) |
| Seat Width: | 0.43m (17in) |
| Seat Height (from the ground): | 0.56m (22in) |
| Battery : | Power Type/ Lithium battery |
| Motor: | 24V 120W |
| Battery Charger: | 29.4V 2A |