

P310 Regal

| | |
|--|---|
|  | <h3>Technical Specs</h3> <p>Max Speed: up to 5mph Range: up to 18 miles Turning Radius: 21" Ground Clearance: 3" Motor: DC 24V, 180W Controller: Dynamic Shark 60A/PG VR2 60A Battery: 12V/U1(34AH) x 2pcs Charger: 5A off-board Gradient: 12° Casters: 8" foam filled tires Drive Wheels: 10" foam filled tires Brakes: intelligent, regenerative, electromagnetic</p> |
| <h3>MSRP</h3> <p>US: \$5999.00</p> | |
| <h3>Medicare Code</h3> <p>K0823 P310-2S C K0822 P310-2S S</p> | |
| <h3>Certification</h3>  | |
| <h3>Features</h3> <ul style="list-style-type: none">• Stylish and compact rear wheel drive model with 10" drive wheels and alloy rims• Extremely low center of gravity• Lower seat height which facilitates transfers and makes it the most stable unit in its class• Sturdy, deluxe captain seat with angle-adjustable back and flip-up armrests• Excellent maneuverability• Programmable controller with a wide range of selections to suit the needs and demands of different users• Non-programmable controller also available• Off-board charger• Semi-reclining seat• Colors Available: Red, Blue, Black• Options: Elevating Legrests, Oxygen Tank Holder, Cane & Crutch Holder, Walker Holder, Rear Foldable Basket | |

Dimensions

P310

| | |
|--|----------|
|  Seat Width | 20 |
|  Seat Depth | 20 |
|  Seat Height | 20-22 |
|  Armrest Height | n/a |
|  Backrest Height | 18 |
|  Min. Footrest Extension | n/a |
|  Max. Footrest Extension | n/a |
|  Overall Length | 39 |
|  Overall Height | 41 |
|  Overall Width | 24 |
|  Folded Width | n/a |
|  Total Weight | 211 lbs. |
| Capacity | 300 lbs. |