

# Junior

## Light Mini Powerbase

The Junior offers stylish design and excellent maneuverability. It is lightweight with swivel seat, flip up full length armrests, headrest and flip-up footplate. Lightweight and easy disassembly, with a removable battery pack for ease of charging.



# Junior

## Features

---

- Lightweight and easy to disassemble for handling storage and transporting
- Ideal as personal secondary portable power wheelchair
- Programmable controller with a wide range of selections to suit the need and demands of different users
- All standard power wheelchair features
- Excellent maneuverability
- Off-board charger

## Specification

Weight Capacity	300lbs
Maximum Speed	4.5 mph
Per-Range Charge	12.5 miles
Batteries	12V/22A x2
Ground Clearance	2"
Turning Radius	23"
Overall Length	35"
Overall Width	33"
Drive Wheels	9" foam filled tires
Casters	7 x 1.75" solid tires
Suspension	No suspension
Seat Width	18"
Seat Depth	18"
Back Height	18"
Warranty	18 months: Motor & Electronics 5 Year/Limited Lifetime Frame
Controller	LinX

## Models

Junior C Rear Wheel Drive	P320-2S C
Junior S Rear Wheel Drive	P320-2S S

## Options

---

- Elevating leg rests
- Oxygen tank holder
- Cane and Crutch holder
- Walker holder
- Foldable rear basket

